


# ENDOLINE® NASO-DRAIN

20

NASOBILIARY CATHETER




The ENDOLINE® NASO-DRAIN nasobiliary catheter is used to aspirate bile or pancreatic juices. This medical device is intended primarily for persons with obstructive jaundice, acute cholangitis or extrahepatic cholestasis.


## ADVANTAGES

- Axial and lateral perforations in the distal part enable efficient drainage of the hepatobiliary tract.
- Easy positioning of the catheter in the bile ducts thanks to the radiopaque distal extremity (with gold ring).

## AVAILABLE FORMS



CURVE 1



CURVE 2



CURVE 6

## AVAILABLE REFERENCE(S)

	Reference	Catheter diameter (Fr)	Useful length (cm)
<b>CURVE 1</b>	AS0530525	5	250
	AS0530725	7	250
	AS0531025	10	250
<b>CURVE 2</b>	AS0550525	5	250
	AS0550725	7	250
<b>CURVE 6</b>	AS0570525	5	250
	AS0570725	7	250

## ABOUT THIS MEDICAL DEVICE

1 <sup>st</sup> release date	1996
Single use	Yes
Class	IIa
Sterile	Yes
Sterilisation method	Ethylene oxide sterilisation
Contains latex	No
Lifetime	5 years
Certification	CE
Packaging	Box of 5 units in individual sterile wrapping
Storage	In its original box in a cool and dry place

PRINCE MEDICAL SAS  
 ZA La Sente Du Moulin  
 64 rue Gutenberg - 60530 Ercuis - France  
[www.prince-medical.com](http://www.prince-medical.com)  
 Tel: +33 (0)3 44 26 54 76  
[prince-medical@omerin.com](mailto:prince-medical@omerin.com)

PRINCE MEDICAL  
 DEVICES

The information provided on this data sheet is indicative and subject to change without notice. PRINCE MEDICAL SAS shall not be held responsible for any incidents resulting from inappropriate use, failure to comply with the rules and regulations in force, or failure to follow the instructions for use provided. Medical devices sold by PRINCE MEDICAL SAS should only be used by a duly qualified practitioner or authorised person.

® Registered trademark of PRINCE MEDICAL SAS. Drawings and photos are not contractually binding. Reproduction is prohibited without the prior consent of PRINCE MEDICAL SAS.