


ENDOLINE® 2V - JECT

28

2 WAYS INJECTION NEEDLE



PLASTIC TIP



METAL TIP

The ENDOLINE® 2V-JECT 2-ways injection needle has been designed specifically and exclusively for haemostasis of gastrointestinal haemorrhages, sclerosis of oesophageal varices and inflating and marking polyps.

The needle is only suitable for use with flexible endoscopes.

The 2nd way can be used to clean cleaning and/or rinse the procedure area.

ADVANTAGES

- Radiopaque distal part visible on the screen, for models with a metal tip.
- Ergonomically shaped handle for easy manipulation.
- Product compatible for use with the main endoscopes on the market.
- Luer-Lock connector for connecting a syringe for product injection.

AVAILABLE HANDLES



1.0
The most widely used handle in the world, used and approved for over 20 years!



2.0
Evolution of the handle 1.0, thicker for a better grip in the hand.

ABOUT THIS MEDICAL DEVICE

1 st release date	1996
Single use	Yes
Class	Ila
Sterile	Yes
Sterilisation method	Ethylene oxide sterilisation
Contains latex	No
Lifetime	5 years
Certification	CE
Packaging	Box of 5 units in individual sterile wrapping
Storage	In its original box in a cool and dry place

PRINCE MEDICAL SAS
ZA La Sente Du Moulin
64 rue Gutenberg - 60530 Ercuis - France
www.prince-medical.com
Tel: +33 (0)3 44 26 54 76
prince-medical@omerin.com



The information provided on this data sheet is indicative and subject to change without notice. PRINCE MEDICAL SAS shall not be held responsible for any incidents resulting from inappropriate use, failure to comply with the rules and regulations in force, or failure to follow the instructions for use provided. Medical devices sold by PRINCE MEDICAL SAS should only be used by a duly qualified practitioner or authorised person.

® Registered trademark of PRINCE MEDICAL SAS. Drawings and photos are not contractually binding. Reproduction is prohibited without the prior consent of PRINCE MEDICAL SAS.

ENDOLINE® 2V - JECT

HANDLE 1.0

> Metal Tip

Reference	Catheter diameter (mm/Fr)	Needle diameter (mm) (G)		Useful length (cm)	Length of the needle outlet (mm)
AS1840718	2,35 / 7	0,7	22	180	5
AS1840723	2,35 / 7	0,7	22	230	5
AS1841018	2,35 / 7	1	19	180	5,5
AS1841023	2,35 / 7	1	19	230	5,5

HANDLE 2.0

> Plastic Tip

AS2830518	1,75/5,3	0,5	25	180	4
AS2830523	1,75/5,3	0,5	25	230	4
AS2830718	2,35/7	0,7	22	180	5
AS2830723	2,35/7	0,7	22	230	5
AS2830728	2,35/7	0,7	22	280	5
AS2831018	2,35/7	1	19	180	5,5
AS2831023	2,35/7	1	19	230	5,5

> Metal Tip

AS2900718	2,35/7	0,7	22	180	5
AS2900723	2,35/7	0,7	22	230	5
AS2900728	2,35/7	0,7	22	280	5
AS2901018	2,35/7	1	19	180	5,5
AS2901023	2,35/7	1	19	230	5,5

PRINCE MEDICAL SAS
 ZA La Sente Du Moulin
 64 rue Gutenberg - 60530 Ercuis - France
www.prince-medical.com
 Tel: +33 (0)3 44 26 54 76
prince-medical@omerin.com