


ENDOLINE® 2V-TRACT

12

2-LUMENS BALLOON CATHETER



QUICK EXCHANGE SYSTEM
Makes it easier to insert catheter
into echo-endoscope
AS2340720

**WELDED
DERIVATIONS**

**BONDED
DERIVATIONS**

TRANSPARENT TUBE
AS4460720

The **ENDOLINE® 2V-TRACT** 2-lumens balloon catheter enables the removal of gallstones and injection of contrast material or any other medicated substance into the bile ducts. This medical device is intended mainly for persons with hepatic colic following the appearance of gallstones.

ADVANTAGES

- Product compatible for use with the main endoscopes on the market.
- Radiopaque.
- Standard connector for connecting a syringe for contrast medium injection.
- Presence of a stainless steel ring and stainless steel ligatures on the balloon reinforcing the radiopacity of the distal part for better balloon localisation.

ADDITIONAL ACCESSORIES

The balloon catheter can be used with the **ENDOLINE® Y-TRACT** Y-clamp fitting. Provides practitioners with 2 injection ports on balloon catheters with a single injection port. It also enables the fixation of the guide wire after catheter placement.

AVAILABLE BALLOONS



LATEX BALLOON CATHETER

PRECAUTIONS FOR USE
**DO NOT USE IN PATIENTS
ALLERGIC TO LATEX**



SILICONE BALLOON CATHETER

for patients
allergic to latex.



ABOUT THIS MEDICAL DEVICE

| | |
|------------------------------|---|
| 1 st release date | 1996 |
| Single use | Yes |
| Class | Ila |
| Sterile | Yes |
| Sterilisation method | Ethylene oxide sterilisation |
| Latex present | Yes except reference AS0920720 |
| Lifetime | 2 years except reference AS0920720 5 years |
| Certification | CE |
| Packaging | Box of 5 units in individual sterile wrapping |
| Storage | In its original box in a cool and dry place |

PRINCE MEDICAL SAS
ZA La Sente Du Moulin
64 rue Gutenberg - 60530 Ercuis - France
www.prince-medical.com
Tel: +33 (0)3 44 26 54 76
prince-medical@omerin.com

PRINCE MEDICAL
DEVICES

The information provided on this data sheet is indicative and subject to change without notice. PRINCE MEDICAL SAS shall not be held responsible for any incidents resulting from inappropriate use, failure to comply with the rules and regulations in force, or failure to follow the instructions for use provided. Medical devices sold by PRINCE MEDICAL SAS should only be used by a duly qualified practitioner or authorised person.

® Registered trademark of PRINCE MEDICAL SAS. Drawings and photos are not contractually binding. Reproduction is prohibited without the prior consent of PRINCE MEDICAL SAS.

ENDOLINE® 2V-TRACT

LATEX
BALLOON

| Reference | Catheter diameter (mm/Fr) | Inflated balloon diameter (mm) | Balloon inflation max gas (mL) | Balloon injection max liquid (mL) | Wire-guide recommended (inch) | Syringe (cc) | Useful length (cm) | Option |
|--------------------------------|------------------------------|--------------------------------------|--------------------------------------|---|-------------------------------------|-----------------|-----------------------|-----------------------|
| > Bonded derivations | | | | | | | | |
| AS0030520 | 1.70 / 5 | 12 | 3.0 | 1.5 | .025 | 3 | 200 | - |
| AS1100720 | 2.35 / 7 1.70 / 5 | 16 | 4.5 | 2.5 | .035 | 5 | 200 | - |
| AS1910720 | 2.35 / 7 | 12 | 3.0 | 1.5 | .035 | 3 | 200 | - |
| AS0030720 | 2.35 / 7 | 16 | 4.5 | 2.5 | .035 | 5 | 200 | - |
| AS4460720 | 2.35 / 7 | 16 | 4.5 | 2.5 | .035 | 5 | 200 | Transparent |
| AS0280720 | 2.35 / 7 | 18 | 5.5 | 3.5 | .035 | 5 | 200 | - |
| > Welded derivations | | | | | | | | |
| AS2340720 | 2.35 / 7 | 16 | 4.5 | 2.5 | .035 | 5 | 200 | Quick exchange system |
| AS2040720 | 2.35 / 7 1.70 / 5 | 16 | 4.5 | 2.5 | .035 | 5 | 200 | Tube cover |
| AS2130720 | 2.35 / 7 | 16 | 4.5 | 2.5 | .035 | 5 | 200 | Tube cover |

SILICONE
BALLOON

| | | | | | | | | |
|--------------------------------|----------|----|-----|-----|------|---|-----|---|
| > Bonded derivations | | | | | | | | |
| AS0920720 | 2.35 / 7 | 16 | 4.5 | 2.5 | .035 | 5 | 200 | - |

PRINCE MEDICAL SAS
 ZA La Sente Du Moulin
 64 rue Gutenberg - 60530 Ercuis - France
www.prince-medical.com
 Tel: +33 (0)3 44 26 54 76
prince-medical@omerin.com