


# CLEANBRUSH® KIT-4, KIT-5

64

CLEANING SWAB KIT – 4 OR 5 PARTS



The **CLEANBRUSH® KIT-4, KIT-5** swab kits are intended for cleaning the working channels, valves and valve cages of flexible endoscopes.

They are suitable for use with the main gastroscopes, colonoscopes, duodenoscopes, enteroscopes, bronchoscopes, cystoscopes, ureteroscopes, nephroscopes, hysteroscopes and naso-laryngoscopes on the market.

## KIT COMPONENTS (SEE THE TABLE FOR THE DETAILED COMPONENTS BY REFERENCE)

- **CLEANBRUSH® ENDO**: cleaning brush for working channel
- **CLEANBRUSH® ENDO** (AN1055002, AN1055023, AN1056523): cleaning brush for the aspiration channel
- **CLEANBRUSH® ENDO** (AN1050123, AN1050138): cleaning brush for air/water channel, and/or Water Jet and/or balloon inflation channel
- **CLEANBRUSH® PISTON**: cleaning brush for valves, valve cages
- **CLEANBRUSH® EXO-B**: brush for external cleaning

## ABOUT THIS MEDICAL DEVICE

1 <sup>st</sup> release date	1999
Single use	Yes
Class	I
Sterile	No (ethylene oxide decontamination)
Contains latex	No
Lifetime	5 years
Certification	CE
Packaging	Box of 25 or 50 units depending on product reference
Storage	In its original box in a cool and dry place

PRINCE MEDICAL SAS  
 ZA La Sente Du Moulin  
 64 rue Gutenberg - 60530 Ercuis - France  
[www.prince-medical.com](http://www.prince-medical.com)  
 Tel: +33 (0)3 44 26 54 76  
[prince-medical@omerin.com](mailto:prince-medical@omerin.com)

**PRINCE MEDICAL**  
 DEVICES

The information provided on this data sheet is indicative and subject to change without notice. PRINCE MEDICAL SAS shall not be held responsible for any incidents resulting from inappropriate use, failure to comply with the rules and regulations in force, or failure to follow the instructions for use provided. Medical devices sold by PRINCE MEDICAL SAS should only be used by a duly qualified practitioner or authorised person.

® Registered trademark of PRINCE MEDICAL SAS. Drawings and photos are not contractually binding. Reproduction is prohibited without the prior consent of PRINCE MEDICAL SAS.

# CLEANBRUSH® KIT-4, KIT-5

**AVAILABLE REFERENCE(S)**

> KIT with 4 parts

Kit reference	REFERENCES FOR CLEANING BRUSHES				OTHER KIT COMPONENTS	Packaging
	For working channel	For aspiration channel	For valves	For air/water channel, water jet, balloon inflation channel	External cleaning brush	
AN1903018	AN1053018 (with 1 brush Ø 3 mm L180 cm)	AN1055023 (with 1 brush Ø 5 mm L230 cm)				50
AN1983018	AN1063018 (with 2 brushes Ø 3 mm L180 cm)			AN1050123 (with 1 brush Ø 1.8 mm L230 cm)		
AN4535018	AN1065018 (with 2 brushes Ø 5 mm L180 cm)	AN1056523 (with 1 brush Ø 6.5 mm L230 cm)				25
AN2543010	AN1053010 (with 1 brush Ø 3 mm L100 cm)	AN1055002 (with 1 brush Ø 5 mm L20 cm)	AN1061110 (with 2 brushes Ø 11 & 5 mm, L15.5 cm)			
AN2553010	AN1063010 (with 2 brushes Ø 3 mm L100 cm)					
AN4505023	AN1055023 (with 1 brush Ø 5 mm L230 cm)					
AN4515023	AN1065023 (with 2 brushes Ø 5 mm L230 cm)					
AN4505018	AN1055018 (with 1 brush Ø 5 mm L180 cm)	AN1056523 (with 1 brush Ø 6.5 mm L230 cm)				
AN4515018	AN1065018 (with 2 brushes Ø 5 mm L180 cm)					
AN2543018	AN1053018 (with 1 brush Ø 3 mm L180 cm)					
AN2545023	AN1055023 (with 1 brush Ø 5 mm L230 cm)					
AN3335023	AN1065023 (with 2 brushes Ø 5 mm L230 cm)					
AN2995023	AN1055023 (with 1 brush Ø 5 mm L230 cm)				AN1050123 (with 1 brush Ø 1.8 mm L230 cm)	
				AN1050138 (with 1 brush Ø 1.8 mm L380 cm)	AN2580000 (U-shaped with soft hair, useful length 2 x 8 cm)	

**AVAILABLE REFERENCE(S)**

> KIT with 5 parts

Kit reference	REFERENCES FOR CLEANING BRUSHES				OTHER KIT COMPONENTS	Packaging
	For working channel	For aspiration channel	For valves	For air/water channel, water jet, balloon inflation channel	External cleaning brush	
AN4545018	AN1055018 (with 1 brush Ø 5 mm L180 cm)					25
AN4555018	AN1065018 (with 2 brushes Ø 5 mm L180 cm)	AN1056523 (with 1 brush Ø 6.5 mm L230 cm)	AN1061110 (with 2 brushes Ø 11 & 5 mm, L15.5 cm)	AN1050123 (with 1 brush Ø 1.8 mm L230 cm)	AN2580000 (U-shaped with soft hair, useful length 2 x 8 cm)	

The information provided on this data sheet is indicative and subject to change without notice. PRINCE MEDICAL SAS shall not be held responsible for any incidents resulting from inappropriate use, failure to comply with the rules and regulations in force, or failure to follow the instructions for use provided. Medical devices sold by PRINCE MEDICAL SAS should only be used by a duly qualified practitioner or authorised person.

® Registered trademark of PRINCE MEDICAL SAS. Drawings and photos are not contractually binding. Reproduction is prohibited without the prior consent of PRINCE MEDICAL SAS.