

PM-LIFE® PRELEVINJEC

30

ASPIRATION/INJECTION DEVICE FOR IUD



The **PM-LIFE® PRELEVINJEC** aspiration/injection device for intrauterine insemination is used to introduce washed and prepared sperm into the uterine cavity without a syringe, as part of artificial intrauterine insemination (IUI) during medically assisted reproduction (MAR).

ADVANTAGES

- Rigid and notched grip tip for a firm grip.
- Atraumatic tip for smooth introduction into the cervical canal.
- Optimal control of injection flow.
- Indicative graduations from 1 to 5 (one graduation corresponding to approximately 100 microlitres).

AVAILABLE REFERENCE(S)

Reference	Inner diameter of the catheter (mm)	Outer diameter of the catheter (mm / Fr)	Inner diameter of the rigid sheath (mm)	Outer diameter of the rigid sheath (mm / Fr)	Length of the catheter (mm)
PS1240623	0.8	2.00 / 6	2.65	3.10 / 9	230

ABOUT THIS MEDICAL DEVICE

1 st release date	2002
Single use	Yes
Class	Is
Sterile	Yes
Sterilisation method	Sterilised via gamma radiation
Contains latex	No
Lifetime	5 years
Certification	CE
Packaging	Box of 25 units in individual sterile wrapping
Storage	In its original box in a cool and dry place

PRINCE MEDICAL SAS
 ZA La Sente Du Moulin
 64 rue Gutenberg - 60530 Ercuis - France
www.prince-medical.com
 Tel: +33 (0)3 44 26 54 76
prince-medical@omerin.com

PRINCE MEDICAL
 DEVICES

The information provided on this data sheet is indicative and subject to change without notice. PRINCE MEDICAL SAS shall not be held responsible for any incidents resulting from inappropriate use, failure to comply with the rules and regulations in force, or failure to follow the instructions for use provided. Medical devices sold by PRINCE MEDICAL SAS should only be used by a duly qualified practitioner or authorised person.

® Registered trademark of PRINCE MEDICAL SAS. Drawings and photos are not contractually binding. Reproduction is prohibited without the prior consent of PRINCE MEDICAL SAS.