

PM-CARE® CHERON

CHERON FORCEPS

27



PM-CARE® CHERON Cheron forceps are used to insert or remove an intrauterine device during gynaecological examinations or procedures. During consultations or gynaecological care, it can be used to grip an aseptic dressing, tampon or compress to apply disinfectant and/or clean the cervix or work area.

PM-CARE® CHERON Cheron forceps can be straight or angled to adapt to the shape of the uterus. The practitioner can choose to use straight or angled forceps according to their personal preferences and ease in the act to be performed.



AVAILABLE REFERENCE(S)

Reference	Forceps model
PS3200025	Straight
PS3210025	Angled



ABOUT THIS MEDICAL DEVICE

1 st release date	2011
Single use	Yes
Class	Is
Sterile	Yes
Sterilisation method	Ethylene oxide sterilisation
Contains latex	No
Lifetime	5 years
Certification	CE
Packaging	Box of 20 units in individual sterile wrapping
Storage	In its original box in a cool and dry place

PRINCE MEDICAL SAS
 ZA La Sente Du Moulin
 64 rue Gutenberg - 60530 Ercuis - France
www.prince-medical.com
 Tel: +33 (0)3 44 26 54 76
prince-medical@omerin.com

The information provided on this data sheet is indicative and subject to change without notice. PRINCE MEDICAL SAS shall not be held responsible for any incidents resulting from inappropriate use, failure to comply with the rules and regulations in force, or failure to follow the instructions for use provided. Medical devices sold by PRINCE MEDICAL SAS should only be used by a duly qualified practitioner or authorised person.

® Registered trademark of PRINCE MEDICAL SAS. Drawings and photos are not contractually binding. Reproduction is prohibited without the prior consent of PRINCE MEDICAL SAS.