

ENDOLINE® HOOD

48

HOOD



The **ENDOLINE® HOOD** fits on the distal end of most endoscopes and protects the cardia, oesophagus and pharynx during endoscopic foreign body removal. It also prevents the foreign body from escaping.

ADVANTAGES

- Bell shape to completely cover the foreign body.
- Easy fixation and compatible with main endoscopes on the market.
- Secured by the attached elastic band.

PRECAUTIONS FOR USE

DO NOT USE IN PATIENTS ALLERGIC TO LATEX

AVAILABLE REFERENCE(S)

Reference	Material	Distal tip diameter (mm)	Proximal tip diameter (mm)	Useful length (mm)
AS0420000	Latex	40	8.3	75

ABOUT THIS MEDICAL DEVICE

1 st release date	1998
Single use	Yes
Class	Is
Sterile	Yes
Sterilisation method	Ethylene oxide sterilisation
Contains latex	Yes
Lifetime	2 years
Certification	CE
Packaging	Box of 10 units in individual sterile wrapping
Storage	In its original box, protected from ultraviolet rays, in a cool and dry place

PRINCE MEDICAL SAS
 ZA La Sente Du Moulin
 64 rue Gutenberg - 60530 Ercuis - France
www.prince-medical.com
 Tel: +33 (0)3 44 26 54 76
prince-medical@omerin.com

PRINCE MEDICAL
 DEVICES

The information provided on this data sheet is indicative and subject to change without notice. PRINCE MEDICAL SAS shall not be held responsible for any incidents resulting from inappropriate use, failure to comply with the rules and regulations in force, or failure to follow the instructions for use provided. Medical devices sold by PRINCE MEDICAL SAS should only be used by a duly qualified practitioner or authorised person.

® Registered trademark of PRINCE MEDICAL SAS. Drawings and photos are not contractually binding. Reproduction is prohibited without the prior consent of PRINCE MEDICAL SAS.