

PM-CARE® 1V-STD

1-LUMEN RIGID HYSTEROGRAPHY AND HYSTEROSONOGRAPHY CATHETER WITH FOAM

19



The PM-CARE® 1V-STD 1-lumen rigid hystero catheter with foam is used to perform hystero and hysterosonography and to ensure a tight seal during the examination. The sliding catheter also enables a deeper penetration into the uterine cavity for injection.

ADVANTAGES

- Easy and simple manipulation.
- Sealing ensured by positioning the foam against the cervix.
- Ease of product injection.
- Atraumatic insertion.

AVAILABLE REFERENCE(S)

Reference	Catheter diameter (mm / Fr)	Foam diameter (mm)	Catheter length (mm)	Sliding catheter
PS1870733	2.45 / 7.5	15	330	-
PS2010733	2.45 / 7.5	25	330	-
PS2020733	1.5 / 4.5	15	340	yes

ABOUT THIS MEDICAL DEVICE

1 st release date	2004 2005 for the reference PS2020733
Single use	Yes
Class	Is
Sterile	Yes
Sterilisation method	Ethylene oxide sterilisation
Contains latex	No
Lifetime	5 years
Certification	CE
Packaging	Box of 5 units in individual sterile wrapping
Storage	In its original box in a cool and dry place



SLIDING

PRINCE MEDICAL SAS
ZA La Sente Du Moulin
64 rue Gutenberg - 60530 Ercuis - France
www.prince-medical.com
Tel: +33 (0)3 44 26 54 76
prince-medical@omerin.com

The information provided on this data sheet is indicative and subject to change without notice. PRINCE MEDICAL SAS shall not be held responsible for any incidents resulting from inappropriate use, failure to comply with the rules and regulations in force, or failure to follow the instructions for use provided. Medical devices sold by PRINCE MEDICAL SAS should only be used by a duly qualified practitioner or authorised person.

® Registered trademark of PRINCE MEDICAL SAS. Drawings and photos are not contractually binding. Reproduction is prohibited without the prior consent of PRINCE MEDICAL SAS.

PRINCE MEDICAL
DEVICES