


ENDOLINE® VIPER

EXTRACTION FORCEPS

41



The ENDOLINE® VIPER double-grasp extraction forceps is used for the extraction of emptied intragastric balloons.

ADVANTAGES

- Three fingers grip for easy manipulation.
- Radiopaque distal part for perfect visibility.
- Product compatible for use with the main endoscopes on the market.
- Reinforced with a sheath spring to support balloon resistance during extraction.

AVAILABLE REFERENCE(S)


Reference	Sheath diameter (mm / Fr)	Sheath length (cm)
AS2290718	2.3/7	180

AVAILABLE AS A KIT WITH THE DRAINAGE NEEDLE**ENDOLINE® PUNC&VIPER / Extraction kit**

Reference: AS2710718

ABOUT THIS MEDICAL DEVICE

1 st release date	2008
Single use	Yes
Class	Is
Sterile	Yes
Sterilisation method	Ethylene oxide sterilisation
Contains latex	No
Lifetime	5 years
Certification	CE
Packaging	Box of 5 units in individual sterile wrapping
Storage	In its original box in a cool and dry place



PRINCE MEDICAL SAS
 ZA La Sente Du Moulin
 64 rue Gutenberg - 60530 Ercuis - France
www.prince-medical.com
 Tel: +33 (0)3 44 26 54 76
prince-medical@omerin.com

The information provided on this data sheet is indicative and subject to change without notice. PRINCE MEDICAL SAS shall not be held responsible for any incidents resulting from inappropriate use, failure to comply with the rules and regulations in force, or failure to follow the instructions for use provided. Medical devices sold by PRINCE MEDICAL SAS should only be used by a duly qualified practitioner or authorised person.

® Registered trademark of PRINCE MEDICAL SAS. Drawings and photos are not contractually binding. Reproduction is prohibited without the prior consent of PRINCE MEDICAL SAS.

PRINCE MEDICAL
 DEVICES